Accutrend® Cholesterol

ES

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EN

EN

the test strips from humidity and direct sunlight. If Accutrend Glucose and Accutrend Cholesterol test strips are used at the same time, please make sure that the container for the glucose test strips is always closed with the white stopper and the container for the cholesterol test strips is closed with the blue stopper. If the container stoppers are mixed up, the specified shelf-life can no longer be guaranteed. Never allow liquids or disinfectants to get into the container. The stopper contains a non-toxic silicate-based desiccant. If the desiccant is inadvertently ingested, please drink plenty of water! All components of the pack can be discarded along with household waste or, when used in a laboratory or a doctor's office, according to the specific guidelines.

Contents of the pack: Tube with 25 test strips, 1 code strip

Additionally required material: One of the above listed meters with instruction manual, lancet device (e.g. Accu-Chek Safe-T-Pro Plus). Performance checks: Accutrend Control CH1 control solution, cellulose swab.

Coding: Whenever a new pack of Accutrend Cholesterol test strips is opened, your meter must be recoded. This procedure is described in the operating instructions of your meter. If the code number on the label of the test strip container and the code number displayed by the meter do not match, the meter will not perform a measurement. Always keep the code strip until all test strips of a pack are used up. However, the code strip should not be stored in the test strip container as this may impair the quality of the test strips.

Checks required everytime you perform a test: Before taking a measurement, the test strip must be checked to see that it is still fit for use. If the window on the back of the test strip has taken on a bluish colour, there is a risk that the cholesterol values displayed are too high. If this occurs, use a new test strip. Every time you perform a test, make sure the test area is completely covered with blood. Follow also the instruction given in the operating instructions of your meter. Occasionally, the round window on the back of the test strip may be mottled; this uneven coloration is sample-related and does not interfere with the measurement of the Accutrend Cholesterol test strip.

Performance check: The performance of the meter and the test strips and their correct handling must be checked regularly using Accutrend Control CH1 control solution. For further details, please refer to the operating instructions for the meter and the pack-age insert of the Accutrend Control CH1 control solution kit. The acceptable control range for a given test strip lot (code number) in combination with the control solution lot is shown on one of the notes supplied with the test strip pack and the control solution pack. The target value corresponds to the mean of the lower and upper value of this control range. Values obtained should fall within the defined limits. Each laboratory must establish corrective measures, if values fall outside the limits. Follow the applicable government regulations and local guidelines for quality control.

Limitations and sources of error: The cholesterol determination may be affected by: intravenous infusion of ascorbic acid (vitamin C); bilirubin values greater than 10 mg/dL (171 µmol/L), e.g. in jaundice; hematocrit values exceeding 55 %; methylaminoantipyrine; gentisic acid.

Measuring range: 150–300 mg/dL (3.88–7.76 mmol/L). Linearity of the method is given in range of 150 mg/dL to 300 mg/dL (3.88–7.76 mmol/L).

Normal values: Cholesterol values less than 200 mg/dL (5.2 mmol/L) are considered normal. In some patient groups, the reference ranges may be different; you should discuss this individually with your doctor.

Notes for doctors and laboratories:

Applications: For the determination of cholesterol in fresh capillary blood in diabetics and non-diabetics; for self-monitoring; for the early detection of a risk of atherosclerosis; for monitoring treatment with lipid-lowering drugs; for screening purposes. If measurements are performed in several persons, as may be the

case in a doctor's office or in screening programmes, the fresh capillary blood may also be applied with heparin-coated capillary pipettes (volume 15 µL to 40 µL).

Test principle: Enzymatic cleavage of the cholesterol ester into fatty acids and cholesterol; oxidation of cholesterol to cholesterol-one with the simultaneous formation of hydrogen peroxide which oxidises an indicator to its blue radical cation.

Components per test: CHE (microorganisms) 1.44 U; CHOD (Nocardia erythropolis) 0.12 U; POD (horseradish) 0.20 U; 3,3',5,5'-tetramethylbenzidine 8.5 µg; non-reactive components 1.4 mg.

Performance data:

Accuracy (method comparison): In method comparisons with capillary blood in hospitals > 95% of all values were found within ±15% of the reference method (CHOD-PAP, Roche Diagnostics). The mean systematic differences were between +2.6% and -3.6%. **Repeatability** (within series imprecision): In investigations with EDTA venous blood in a concentration range between 167 mg/dL (4.32 mmol/L) and 291 mg/dL (6.75 mmol/L) typical measurement series showed a CV (coefficient of variation) of 0.8 to 3.7%. **Reproducibility** (day-to-day imprecision): The day-to-day imprecision was determined in test series with control solution. At a concentration of 191 mg/dL (approx. 5 mmol/L) the CV was typically between 1.1% and 3.8%, at a concentration of approx. 270 mg/dL (ca. 7 mmol/L) between 2.4% and 5.0%. **Detection limit** (lowest value detected): 150 mg/dL (3.88 mmol/L).

Last revision: 2010-03

Presentation: 25 test strips and 1 code strip **REF** 11418262171

For literature references and an explanation of the symbols used please refer to the end of this package insert.

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CE 0123

ACCUTREND, ACCU-CHEK and COBAS are trademarks of Roche.

© 2010 Roche Diagnostics GmbH.

 Roche Diagnostics GmbH
Sandhofer Strasse 116
D-68305 Mannheim
www.roche.com
www.diavant.com
Made in Germany

05956307001(01)
R1 (black)